

Product Datasheet

Human C-type Lectin Domain Family 4 Member C ELISA Kit (orb1213277)

Description	Human C-type Lectin Domain Family 4 Member C, CLEC4C ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	CLEC4C
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,038
Uniprot ID	Q8WTT0
Solubility (25°C)	48ng/ml
Sensitivity	0.054ng/ml
Entrez	170482
Expiration Date	Please enquire.