

Product Datasheet

Human Complement Component C6 ELISA Kit (orb1213288)

Description	Human Complement Component C6, C6 ELISA Kit
Reactivity	Human
Range	2-700ug/ml
Tested Applications	ELISA
Target	C6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	104,786
Uniprot ID	P13671
Solubility (25°C)	9600µg/ml
Sensitivity	1.17µg/ml
Entrez	729
Expiration Date	Please enquire.