

Product Datasheet

Human Complement C1Q Subcomponent Subunit C ELISA Kit (orb1213291)

Description	Human Complement C1Q Subcomponent Subunit C, C1QC ELISA Kit
Reactivity	Human
Range	5-1000ng/mL
Tested Applications	ELISA
Target	C1QC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,774
Uniprot ID	P02747
Solubility (25°C)	1280ng/mL
Sensitivity	2.49ng/mL
Entrez	714
Expiration Date	Please enquire.