

Product Datasheet

Human Complement C1Q Subcomponent Subunit A ELISA Kit (orb1213292)

Description	Human Complement C1Q Subcomponent Subunit A, C1QA ELISA Kit
Reactivity	Human
Range	5-1000ng/mL
Tested Applications	ELISA
Target	C1QA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,017
Uniprot ID	P02745
Solubility (25°C)	1280ng/mL
Sensitivity	2.45ng/mL
Entrez	712
Expiration Date	Please enquire.