

Product Datasheet

Human Collectin-12 ELISA Kit (orb1213293)

Description	Human Collectin-12, COLEC12 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	COLEC12
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	81,515
Uniprot ID	Q5KU26
Solubility (25°C)	48ng/ml
Sensitivity	0.043ng/ml
Entrez	81035
Expiration Date	Please enquire.