

Product Datasheet

Human Cold-inducible RNA-binding Protein ELISA Kit (orb1213294)

Description	Human Cold-inducible RNA-binding Protein, CIRBP ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	CIRBP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,648
Uniprot ID	Q14011
Solubility (25°C)	1600ng/L
Sensitivity	4.02ng/L
Entrez	1153
Expiration Date	Please enquire.