

Product Datasheet

Human Chordin-like Protein 1 ELISA Kit (orb1213303)

Description	Human Chordin-like Protein 1, CHRD L1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CHRD L1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,026
Uniprot ID	Q9BU40
Solubility (25°C)	24ng/ml
Sensitivity	0.027ng/ml
Entrez	91851
Expiration Date	Please enquire.