

Product Datasheet

Human Chondroitin Sulfate Proteoglycan 4 ELISA Kit (orb1213305)

Description	Human Chondroitin Sulfate Proteoglycan 4, CSPG4 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CSPG4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	250,537
Uniprot ID	Q6UVK1
Solubility (25°C)	3200ng/L
Sensitivity	6.91ng/L
Entrez	1464
Expiration Date	Please enquire.