

Product Datasheet

Human Cholesteryl Ester Transfer Protein ELISA Kit (orb1213308)

Description	Human Cholesteryl Ester Transfer Protein, CETP ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	CETP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,756
Uniprot ID	P11597
Solubility (25°C)	96ng/mL
Sensitivity	0.14ng/mL
Entrez	1071
Expiration Date	Please enquire.