

Product Datasheet

Human Chloride Intracellular Channel Protein 4 ELISA Kit (orb1213309)

Description	Human Chloride Intracellular Channel Protein 4, CLIC4 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CLIC4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,772
Uniprot ID	Q9Y696
Solubility (25°C)	3200ng/L
Sensitivity	7.15ng/L
Entrez	25932
Expiration Date	Please enquire.