

Product Datasheet

Human Centromere Protein K ELISA Kit (orb1213315)

Description	Human Centromere Protein K, CENPK ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CENPK
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	31,655
Uniprot ID	Q9BS16
Solubility (25°C)	3200ng/L
Sensitivity	6.95ng/L
Entrez	64105
Expiration Date	Please enquire.