

Product Datasheet

Human Cd99 Antigen-like Protein 2 ELISA Kit (orb1213319)

Description	Human Cd99 Antigen-like Protein 2, CD99L2 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CD99L2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,986
Uniprot ID	Q8TCZ2
Solubility (25°C)	24ng/ml
Sensitivity	0.028ng/ml
Entrez	83692
Expiration Date	Please enquire.