

Product Datasheet

Human Cca Trna Nucleotidyltransferase 1 ELISA Kit (orb1213323)

Description	Human Cca Trna Nucleotidyltransferase 1, Mitochondrial, TRNT1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	TRNT1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,128
Uniprot ID	Q96Q11
Solubility (25°C)	3200ng/L
Sensitivity	7.49ng/L
Entrez	51095
Expiration Date	Please enquire.