

Product Datasheet

Human Cytidine Deaminase ELISA Kit (orb1213327)

Description	Human Cytidine Deaminase, CDA ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	CDA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,185
Uniprot ID	P32320
Solubility (25°C)	2400ng/L
Sensitivity	5.43ng/L
Entrez	978
Expiration Date	Please enquire.