

Product Datasheet

Human Argininosuccinate Synthase ELISA Kit (orb1213335)

Description	Human Argininosuccinate Synthase, ASS1 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	ASS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,530
Uniprot ID	P00966
Solubility (25°C)	960ng/L
Sensitivity	1.21ng/L
Entrez	445
Expiration Date	Please enquire.