

Product Datasheet

Human Apolipoprotein D ELISA Kit (orb1213338)

Description	Human Apolipoprotein D, APOD ELISA Kit
Reactivity	Human
Range	3-900ng/mL
Tested Applications	ELISA
Target	APOD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,276
Uniprot ID	P05090
Solubility (25°C)	960ng/mL
Sensitivity	1.33ng/mL
Entrez	347
Expiration Date	Please enquire.