

Product Datasheet

Human Angiogenic Factor with G Patch and Fha Domains 1 ELISA Kit (orb1213345)

Description	Human Angiogenic Factor with G Patch and Fha Domains 1, AGGF1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	AGGF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	80,977
Uniprot ID	Q8N302
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	55109
Expiration Date	Please enquire.