

Product Datasheet

Human An1-type Zinc Finger Protein 6 ELISA Kit (orb1213347)

Description	Human An1-type Zinc Finger Protein 6, ZFAND6 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	ZFAND6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,555
Uniprot ID	Q6FIF0
Solubility (25°C)	3200ng/L
Sensitivity	7.65ng/L
Entrez	54469
Expiration Date	Please enquire.