

Product Datasheet

Human Amyloid-like Protein 1 ELISA Kit (orb1213349)

Description	Human Amyloid-like Protein 1, APLP1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	APLP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	72,176
Uniprot ID	P51693
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	333
Expiration Date	Please enquire.