

Product Datasheet

Human Alpha-N-acetylgalactosaminidase ELISA Kit (orb1213352)

Description	Human Alpha-N-acetylgalactosaminidase, NAGA ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	NAGA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,565
Uniprot ID	P17050
Solubility (25°C)	96ng/mL
Sensitivity	0.14ng/mL
Entrez	4668
Expiration Date	Please enquire.