

Product Datasheet

Human Alpha-crystallin B Chain ELISA Kit (orb1213354)

Description	Human Alpha-crystallin B Chain, CRYAB ELISA Kit
Reactivity	Human
Range	7-1500.ng/L
Tested Applications	ELISA
Target	CRYAB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,159
Uniprot ID	P02511
Solubility (25°C)	1600ng/L
Sensitivity	3.36ng/L
Entrez	1410
Expiration Date	Please enquire.