

Product Datasheet

Human Alpha-2-Macroglobulin-like Protein 1 ELISA Kit (orb1213357)

Description	Human Alpha-2-Macroglobulin-like Protein 1, A2ML1 ELISA Kit
Reactivity	Human
Range	5-1000ng/mL
Tested Applications	ELISA
Target	A2ML1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	161,107
Uniprot ID	A8K2U0
Solubility (25°C)	1280ng/mL
Sensitivity	2.29ng/mL
Entrez	144568
Expiration Date	Please enquire.