

Product Datasheet

Human Adp-Ribosylation Factor-like Protein 8B ELISA Kit (orb1213361)

Description	Human Adp-Ribosylation Factor-like Protein 8B, ARL8B ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	ARL8B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,539
Uniprot ID	Q9NVJ2
Solubility (25°C)	2400ng/L
Sensitivity	5.14ng/L
Entrez	55207
Expiration Date	Please enquire.