

Product Datasheet

Human Adenosylhomocysteinase ELISA Kit (orb1213363)

Description	Human Adenosylhomocysteinase, AHCY ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	AHCY
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,716
Uniprot ID	P23526
Solubility (25°C)	4800ng/L
Sensitivity	9.97ng/L
Entrez	191
Expiration Date	Please enquire.