

Product Datasheet

Human Actin ELISA Kit (orb1213368)

Description	Human Actin, Alpha Cardiac Muscle 1, ACTC1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	ACTC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,019
Uniprot ID	P68032
Solubility (25°C)	3200ng/L
Sensitivity	7.45ng/L
Entrez	70
Expiration Date	Please enquire.