

Product Datasheet

Human 3-Hydroxyacyl-CoA Dehydrogenase Type-2 ELISA Kit (orb1213379)

Description	Human 3-Hydroxyacyl-CoA Dehydrogenase Type-2, HSD17B10 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	HSD17B10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,923
Uniprot ID	Q99714
Solubility (25°C)	3200ng/L
Sensitivity	7.19ng/L
Entrez	3028
Expiration Date	Please enquire.