

Product Datasheet

Human Low Affinity Immunoglobulin Gamma Fc Region Receptor 3-A ELISA Kit (orb1213391)

Description	Human Low Affinity Immunoglobulin Gamma Fc Region Receptor 3-A, FCGR3A ELISA Kit
Reactivity	Human
Range	5-1000ng/ml
Tested Applications	ELISA
Target	FCGR3A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,089
Uniprot ID	P08637
Solubility (25°C)	1280ng/ml
Sensitivity	2.86ng/ml
Entrez	2214
Expiration Date	Please enquire.