

Product Datasheet

Human DNA Damage-inducible Transcript 3 Protein ELISA Kit (orb1213396)

Description	Human DNA Damage-inducible Transcript 3 Protein, DDIT3 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	DDIT3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,175
Uniprot ID	P35638
Solubility (25°C)	64ng/ml
Sensitivity	0.097ng/ml
Entrez	1649
Expiration Date	Please enquire.