

Product Datasheet

Human Zinc Transporter Zip14 ELISA Kit (orb1213410)

Description	Human Zinc Transporter Zip14, ZIP14 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	SLC39A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,250
Uniprot ID	Q9NY26
Solubility (25°C)	6400ng/L
Sensitivity	9.88ng/L
Entrez	27173
Expiration Date	Please enquire.