

Product Datasheet

Human G Protein Coupled Estrogen Receptor 1 ELISA Kit (orb1213412)

Description	Human G Protein Coupled Estrogen Receptor 1, GPER ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	GPER1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,248
Uniprot ID	Q99527
Solubility (25°C)	64ng/ml
Sensitivity	0.099ng/ml
Entrez	2852
Expiration Date	Please enquire.