

Product Datasheet

Human Laminin Subunit Alpha-1 ELISA Kit (orb1213421)

Description	Human Laminin Subunit Alpha-1, LAMA1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	LAMA1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	337,084
Uniprot ID	P25391
Solubility (25°C)	24ng/ml
Sensitivity	0.025ng/ml
Entrez	284217
Expiration Date	Please enquire.