

Product Datasheet

Human Interleukin-1 Family Member 10 ELISA Kit (orb1213442)

Description	Human Interleukin-1 Family Member 10, IL1F10 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	IL1F10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,943
Uniprot ID	Q8WWZ1
Solubility (25°C)	4800ng/L
Sensitivity	11.15ng/L
Entrez	84639
Expiration Date	Please enquire.