

Product Datasheet

Human Tropomyosin 3 ELISA Kit (orb1213446)

Description	Human Tropomyosin 3, TPM3 ELISA Kit
Reactivity	Human
Range	0.1-40.ng/ml
Tested Applications	ELISA
Target	TMOD3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,595
Uniprot ID	Q9NYL9
Solubility (25°C)	48ng/ml
Sensitivity	0.054ng/ml
Entrez	29766
Expiration Date	Please enquire.