

Product Datasheet

Human Cofilin-1 ELISA Kit (orb1213456)

Description	Human Cofilin-1, CFL1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CFL1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,502
Uniprot ID	P23528
Solubility (25°C)	3200ng/L
Sensitivity	7.96ng/L
Entrez	1072
Expiration Date	Please enquire.