

Product Datasheet

Human Troponin I ELISA Kit (orb1213466)

Description	Human Troponin I, Slow Skeletal Muscle, SLOW SKELETAL MUSCLE
Reactivity	Human
Range	5-50ng/ml
Tested Applications	ELISA
Target	TNNI1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,692
Uniprot ID	P19237
Solubility (25°C)	64ng/ml
Sensitivity	2.43ng/ml
Entrez	7135
Expiration Date	Please enquire.