

## Product Datasheet

# Human Methionine Adenosyltransferase 2 Subunit Beta ELISA Kit (orb1213470)

|                            |                                                                      |
|----------------------------|----------------------------------------------------------------------|
| <b>Description</b>         | Human Methionine Adenosyltransferase 2 Subunit Beta, MAT2B ELISA Kit |
| <b>Reactivity</b>          | Human                                                                |
| <b>Range</b>               | 10-2000ng/L                                                          |
| <b>Tested Applications</b> | ELISA                                                                |
| <b>Target</b>              | MAT2B                                                                |
| <b>Storage</b>             | 2-8°C                                                                |
| <b>Note</b>                | For research use only                                                |
| <b>Application notes</b>   | Sandwich                                                             |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                               |
| <b>MW</b>                  | 37,552                                                               |
| <b>Uniprot ID</b>          | <b>Q9NZL9</b>                                                        |
| <b>Solubility (25°C)</b>   | 2400ng/L                                                             |
| <b>Sensitivity</b>         | 5.28ng/L                                                             |
| <b>Entrez</b>              | <b>27430</b>                                                         |
| <b>Expiration Date</b>     | Please enquire.                                                      |