

Product Datasheet

Human Multimerin-2 ELISA Kit (orb1213471)

Description	Human Multimerin-2, MMRN2 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	MMRN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	104,409
Uniprot ID	Q9H8L6
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	79812
Expiration Date	Please enquire.