

Product Datasheet

Human Gamma-Glutamyl Hydrolase ELISA Kit (orb1213478)

Description	Human Gamma-Glutamyl Hydrolase, GGH ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	GGH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,964
Uniprot ID	Q92820
Solubility (25°C)	64ng/ml
Sensitivity	0.094ng/ml
Entrez	8836
Expiration Date	Please enquire.