

Product Datasheet

Human Branched-chain-amino-acid Aminotransferase ELISA Kit (orb1213479)

Description	Human Branched-chain-amino-acid Aminotransferase, Mitochondrial, BCAT2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	BCAT2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,288
Uniprot ID	015382
Solubility (25°C)	3200ng/L
Sensitivity	8.01ng/L
Entrez	587
Expiration Date	Please enquire.