

Product Datasheet

Human Aspartyl/Asparaginyl Beta-hydroxylase ELISA Kit (orb1213480)

Description	Human Aspartyl/Asparaginyl Beta-hydroxylase, ASPH ELISA Kit
Reactivity	Human
Range	2-600ng/L
Tested Applications	ELISA
Target	ASPH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	85,863
Uniprot ID	Q12797
Solubility (25°C)	640ng/L
Sensitivity	1.24ng/L
Entrez	444
Expiration Date	Please enquire.