

Product Datasheet

Human Period Circadian Protein Homolog 2 ELISA Kit (orb1213481)

Description	Human Period Circadian Protein Homolog 2, PER2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	PER2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	136,579
Uniprot ID	O15055
Solubility (25°C)	3200ng/L
Sensitivity	7.37ng/L
Entrez	8864
Expiration Date	Please enquire.