

Product Datasheet

Human Moesin ELISA Kit (orb1213495)

Description	Human Moesin, MSN ELISA Kit
Reactivity	Human
Range	200-4500ng/L
Tested Applications	ELISA
Target	MSN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	67,820
Uniprot ID	P26038
Solubility (25°C)	4800ng/L
Sensitivity	9.47ng/L
Entrez	4478
Expiration Date	Please enquire.