

Product Datasheet

Human Thymosin Beta 4 ELISA Kit (orb1213498)

Description	Human Thymosin Beta 4, TMSB4 ELISA Kit
Reactivity	Human
Range	30-9000ng/ml
Tested Applications	ELISA
Target	TMSB4X
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	5,053
Uniprot ID	P62328
Solubility (25°C)	9600ng/ml
Sensitivity	15.79ng/ml
Entrez	7114
Expiration Date	Please enquire.