

Product Datasheet

Human Plastin-2 ELISA Kit (orb1213502)

Description	Human Plastin-2, LCP1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	LCP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	70,288
Uniprot ID	P13796
Solubility (25°C)	4800ng/L
Sensitivity	9.62ng/L
Entrez	3936
Expiration Date	Please enquire.