

Product Datasheet

Human Calcium Transporting Atpase Type 2C Member 1 ELISA Kit (orb1213507)

Description	Human Calcium Transporting Atpase Type 2C Member 1, ATP2C1ELISA Kit
Reactivity	Human
Range	1-16ng/mL
Tested Applications	ELISA
Target	ATP2C1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	100,577
Uniprot ID	P98194
Solubility (25°C)	32ng/mL
Sensitivity	0.54ng/mL
Entrez	27032
Expiration Date	Please enquire.