

Product Datasheet

Human C-C Motif Chemokine 25 ELISA Kit (orb1213508)

Description	Human C-C Motif Chemokine 25, CCL25 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	CCL25
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,609
Uniprot ID	O15444
Solubility (25°C)	1600ng/L
Sensitivity	3.43ng/L
Entrez	6370
Expiration Date	Please enquire.