

Product Datasheet

Human Osteocrin ELISA Kit (orb1213509)

Description	Human Osteocrin, OSTN ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	OSTN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,722
Uniprot ID	P61366
Solubility (25°C)	24ng/ml
Sensitivity	0.023ng/ml
Entrez	344901
Expiration Date	Please enquire.