

Product Datasheet

Human Cytoplasmic Aconitate Hydratase ELISA Kit (orb1213519)

Description	Human Cytoplasmic Aconitate Hydratase, ACO1 ELISA Kit
Reactivity	Human
Range	2-32ng/mL
Tested Applications	ELISA
Target	ACO1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	98,399
Uniprot ID	P21399
Solubility (25°C)	64ng/mL
Sensitivity	1.03ng/mL
Entrez	48
Expiration Date	Please enquire.