

Product Datasheet

Human Pannexin-1 ELISA Kit (orb1213526)

Description	Human Pannexin-1, PANX1 ELISA Kit
Reactivity	Human
Range	0.5-8ng/mL
Tested Applications	ELISA
Target	PANX1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	48,050
Uniprot ID	Q96RD7
Solubility (25°C)	16ng/mL
Sensitivity	0.021ng/mL
Entrez	24145
Expiration Date	Please enquire.