

Product Datasheet

Human Properdin ELISA Kit (orb1213527)

Description	Human Properdin, CFP ELISA Kit
Reactivity	Human
Range	7-1500ng/ml
Tested Applications	ELISA
Target	CFP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,276
Uniprot ID	P27918
Solubility (25°C)	1600ng/ml
Sensitivity	3.82ng/ml
Entrez	5199
Expiration Date	Please enquire.